

Genomic DNA Extraction, Amplification and T7E1 Assay

T7 endonuclease cleaves double-stranded DNA at positions of mismatches. NHEJ repair of Cas9-induced breaks will leave a variety of different mutations, and there will almost always be some wild type sequence remaining. Thus, when you amplify the target region, denature and renature the products, there will be mismatches at the target site, if cleavage was effective.

Genomic DNA extraction:

Trypsinize the cells, and put half of the volume in a 1.5 ml microfuge tube and spin briefly at 4°.

Remove supernate and wash cells with 1 ml PBS.

Add 100 µl of Quick Extract

Vortex 15 sec

Heat at 65° for 6 min, then 2 min at 95° (programmed for you)

Measure DNA concentration with Nanodrop spectrophotometer

Dilute to 20 ng/µl with more Quick Extract

PCR amplification of target:

Mix:		Final conc.
29.0 µl	H ₂ O	
10.0 µl	5x Q5 Buffer	1x
1.0 µl	10 mM dNTPs	200 µM
2.5 µl	Fwd primer (10 µM)	0.5 µM
2.5 µl	Rev primer (10 µM)	0.5 µM
5.0 µl	Genomic DNA	100 ng
0.5 µl	Q5 DNA polymerase (2u/µl)	1 unit
50.0 µl		

PCR: 98° 30 sec
98° 10 sec
66° 30 sec } 30x
72° 20 sec
72° 2 min
4° hold

Run 3 µl of the PCR product on a 2% agarose gel.

Purification of PCR product:

Add 250 µl of Zymo DNA binding buffer to the PCR product

Load onto Zymo Clean and Concentrate spin column

Spin 30 sec in microfuge; discard flow through

Add 200 µl DNA Wash buffer; spin 30 sec; discard flow through

Repeat wash step

Place spin column in a fresh 1.5 ml tube

Add 20 µl of EB (elution buffer); spin 1 min in microfuge

Measure DNA concentration with Nanodrop spectrophotometer

T7 endonuclease assay:

Mix: 10 μ l PCR product
2 μ l 10x NEB Buffer 2
8 μ l H₂O

Anneal: 95° 5 min
95→85° at 2°/sec
85→25° at 0.1°/sec
4° hold

Remove 10 μ l to serve as undigested sample.

To the remainder, add 1 μ l T7 Endonuclease; incubate 15 min, 37°.
Run the digested and undigested samples in adjacent lanes on a 2% agarose gel.

Dana Carroll, June 24, 2016